

3-PLY FACE MASK



DESCRIPTION

- Pleated 3 layer design, the filter efficiency is higher than 98%
- Soft lined, water resistant, non-woven material
- Latex and fiberglass free
- Adjustable nose piece and ear-loop design for convenience and comfort
- Certified at both Levels 2 and 3

3-PLY FACE MASK

OTHER INFORMATION

Material	Polypropylene Non-woven fabric
Filter Rate	BFE ≥ 98%
Use Range	Hospitals/medical (Personal protection)
Disposable	Yes
Type	Non-sterility (Medical)
Style	Ear-loop
Colour	Blue

